

**Product Name :**  
Cryo Vial

**Product Code :**  
CHEMISTRYSUK-AS03LAB010R



**Description :**

Cryo Vial

**Technical Specification :**

The cap of these vials is provided with a seal ring inside to make it completely leakproof.

Cryo vials are normally used for storage of Biological material, human & animal cells at temperature as low as  $-190^{\circ}\text{C}$ .

These sterilized tubes are provided with white, prominent graduations and a long marking area for sample identification.

The longitudinal grooves in the base give a non twisting feature to these vials when placed in a cryo rack.

**Capacity:** 4.5ml, 1.8ml, 1.0ml



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265