

Product Name :
Condensers, Liebig, Borosilicate Glass

Product Code :
CHEMISTRYSUK-AS03LAB008S



Description :

Condensers, Liebig, Borosilicate Glass

Technical Specification :

These Condensers, Liebig, Borosilicate Glass available in different colour can easily fit the cap of Condensers, Liebig, Borosilicate Glass for colour identification of different samples.

Type Capacity

0.01 250mm

0.02 300mm

0.03 400mm

0.04 500mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265