

Product Name :
Pneumatic Troughs, Plastic

Product Code :
CHEMISTRYSUK-AS03LAB008U



Description :

Pneumatic Troughs, Plastic

Technical Specification :

These Pneumatic Troughs, Plastic available in different colour can easily fit the cap of Pneumatic Troughs, Plastic for colour identification of different samples.

Type Capacity

0.01 Dia-180 x h-60mm

0.02 Dia-200 x h-100mm

0.03 Dia-250 x h-100mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265