

Product Name :
Pipette, Volumetric

Product Code :
CHEMISTRYSUK-AS03LAB009U



Description :

Pipette, Volumetric

Technical Specification :

These Pipette, Volumetric available in different colour can easily fit the cap of Pipette, Volumetric for colour identification of different samples.

Type Capacity ml

0.01 1
0.02 2
0.03 5
0.04 10
0.05 20
0.06 25
0.07 50
0.08 100



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265