

Product Name :
Test tube stand, Wire pattern

Product Code :
CHEMISTRYSUK-AS03LAB007W



Description :

Test tube stand, Wire pattern

Technical Specification :

These Test tube stand, Wire pattern available in different colour can easily fit the cap of Test tube stand, Wire pattern for colour identification of different samples.

Capacity

~~031m~~
x 90 Tubes

~~062m~~
x 60 Tubes

~~008m~~
x 40 Tubes

~~054m~~
x 24 Tubes

~~005m~~
x 21 Tubes



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265