

**Product Name :**  
Thermometer Clamp

**Product Code :**  
CHEMISTRYSUK-AS03LAB006A



**Description :**

Thermometer Clamp

**Technical Specification :**

Made of a polished aluminum muff with brass clamping jaws.

Excellent clamp for supporting thermometers.

Jaws can be rotated and locked at any angle.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

