

Product Name :
Boss Head, With Detachable Attachment

Product Code :
CHEMISTRYUK-AS03LAB007A



Description :

Boss Head, With Detachable Attachment

Technical Specification :

A unique design with phosphorus bronze clamps for holding objects vertically or horizontally, gives it more flexibility.

Available with Thermometer Hanger.

Multipurpose Aluminum section boss head with brass headed nickel plated steel thumb screw and plated wing shaped nut for fixing the different attachments like thermometer holders, suspension hangers etc.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265