

**Product Name :**  
Universal Clamp

**Product Code :**  
CHEMISTRYUK-AS03LAB008A



**Description :**

Universal Clamp

**Technical Specification :**

Jaws can clamp cylindrical objects from 9.5 to 35mm diameter.

Universal Clamp comprises of Prong Clamps 8mm in diameter and 150mm long.

Jaws can be rotated through and locked in any position.

Good reach during practical.

Ensures good grip.

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



**Astra Scientific**