


**Product Name :**  
Test Tube Rack With Drying Pegs

**Product Code :**  
CHEMISTRYSUK-AS03LAB005B



**Description :**

Test Tube Rack With Drying Pegs

**Technical Specification :**

Vertical 12 pegs are provided for quick drying of the test tubes.

Giving you more room on a single platform.

A newly designed test tube rack with pegs made of beach wood, well polished, with two row of holes.

We can also customize according to individuals preference.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265